

# Test sample of cellulose nanofibers prepared by ACC

We can provide the test samples of cellulose nanofiber prepared by ACC method.

※ACC method: the aqueous counter collision method

## Types of products

Description	Concentration	Content
<b>CNF-1</b>	ca. 1 wt.%	10 kg
<b>CNF-10</b>	ca. 10 wt.%	200 g

## Raw material pulp type

Cat. no	
<b>BB</b>	Bamboo bleached kraft pulp
<b>LB</b>	Hardwood bleached kraft pulp
<b>NB</b>	Softwood bleached kraft pulp

## Preparation condition

Cat.no	Degree of nano-pulverization	Degree of variation in fiber width
<b>A</b>	Low	Wide
<b>B</b>	Middle	Middle
<b>C</b>	High	Narrow

Chuetsu Pulp & Paper Co., Ltd.

Head office in Takaoka

Development Headquarters

282, Yonejima, Takaoka, Toyama 933-8533

⇒ <https://www.cpc-cenf.com/form.html>

